

SPECIFICATIONS

ROMED POLYETHYLENE APRONS

REF: 1/60, 21/60, 21/60B



Medical device conform MDR (EU) 2017/745

REFERENCE

1/60	White, packed per 100 pieces in a polybag
21/60	White, packed per piece in a polybag
21/60B	Blue, packed per piece in a polybag

MATERIAL

100% low density polyethylene (LDPE) film, embossed

DIMENSIONS

Contact thickness	0.025 mm
Length	125 cm ($\pm 5\%$)
Width	80 cm ($\pm 5\%$)
Weight	16 g ($\pm 5\%$)

These products are medical grade personal protection disposable devices and therefore are also suitable for use in the food industry.

PACKING

21/60 and 21/60B packed per piece in a bag, 100 bags in another bag and 10 bags in an export carton.
1/60 packed per 100 pieces in a polybag and 10 polybags in an export carton.