

Technical Datasheet

Date: 22.06.20

Page 1 of 1

EN ISO 13485 certified		
Article - No:	04-093-0396	Revision: 1
Remarks	fitting plate: 04-072-0120 and 04-072-0140	single use only

Specification

Product	Silicone sealing mat, for 384 deep well plates, square and pre-slit wells, natural, np pcr ready autoclavable up to 121°C
Dimensions	ANSI/SLAS 1
Colour	natural
Material Characteristics	silicone (SI)
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 10
	pcs per carton: 50
	pcs per pallet: 5.850
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	