Clin-Tech Ltd.

61027x
Perls' Solution 2
Revision 05 – Date 09/2015
According to regulation (EC) # 1907/2006

Safety Data Sheet

1. IDENTIFICATION OF THE SUBSTANCE / PREPARATION AND COMPANY

PRODUCT NAME PERLS' SOLUTION 2

PRODUCT CODE 61027x

SUPPLIER Clin-Tech Limited Telephone +44 (0)1483 301902

Unit G, Perram Works Fax +44 (0)1483 301907

Merrow Lane

Guildford, Surrey Email info@clin-tech.co.uk GU4 7BN Website www.clin-tech.co.uk

UK

2. HAZARDS IDENTIFICATION

<u>CLASSIFICATION ACCORDING TO REGULATION (EC) NO 1272/2008 [EU-GHS/CLP]</u> Not a dangerous substance according to GHS.

<u>LABELLING ACCORDING TO REGULATION (EC) NO 1272/2008 [EU-GHS/CLP]</u>
No symbols required

Hazard statement(s)

None required

Precautionary statement(s)

None required

3. COMPOSITION / INFORMATION ON INGREDIENTS

Name	EC#	CAS#	Content	Classification
HYDROCHLORIC ACID	231-595-7	7647-01-0	1-5%	Skin Corr. 1B; STOT SE 3 H314 H335
				11314 11333

4. FIRST-AID MEASURES

INHALATION

Move the exposed person to fresh air. Perform artificial respiration if breathing has stopped.

INGESTION

Rinse mouth with water. Never make an unconscious person vomit or drink fluids.

SKIN CONTACT

Remove contaminated clothes and rinse skin thoroughly with water.

EYE CONTACT

Make sure to remove any contact lenses from the eyes before rinsing. Promptly wash eyes with plenty of water while lifting the eyelids.

61027x Perls' Solution 2 Revision 05 – Date 09/2015

5. FIRE-FIGHTING MEASURES

EXTINGUISHING MEDIA

Use: Water spray, Alcohol-resistant foam, Dry chemicals, or Carbon dioxide

PROTECTIVE MEASURES IN FIRE

Self-contained breathing apparatus if necessary.

6. ACCIDENTAL RELEASE MEASURES

PERSONAL PRECAUTIONS

Wear protective clothing as described in Section 8 of this safety data sheet. Evacuate personnel to a safe area. Avoid aerosol formation. Avoid breathing aerosols. Ensure adequate ventilation.

ENVIRONMENTAL PRECAUTIONS

Do not allow spillage to enter watercourses, drains or sewers. Stop further spillage if possible without risk.

SPILL CLEAN-UP METHODS

Flush area with water – do not allow water to enter drains. Absorb with disposable towels and hold for hazardous disposal – see section 13.

7. HANDLING & STORAGE

USAGE PRECAUTIONS

Avoid spilling, skin and eye contact. Avoid aerosol formation. Provide adequate exhaust function where aerosol is likely.

STORAGE PRECAUTIONS

Store in tightly closed original container in a cool, dry well-ventilated place.

STORAGE CLASS

Chemical storage.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

ENGINEERING MEASURES

Provide adequate general and local exhaust ventilation.

PROTECTIVE EQUIPMENT







Clin-Tech Ltd.

61027x Perls' Solution 2 Revision 05 – Date 09/2015

HAND PROTECTION

Use protective gloves. The most suitable glove should be chosen in consultation with the gloves supplier, who can inform about the breakthrough time of the glove material.

EYE PROTECTION

Wear approved, tight-fitting safety glasses where splashing is probable.

RESPIRATORY EQUIPMENT

Self-contained breathing apparatus where aerosols are likely.

OTHER PROTECTION

Provide eyewash station.

HYGIENE MEASURES

DO NOT SMOKE IN WORK AREA! Wash at the end of each work shift and before eating, smoking and using the toilet. Promptly remove any clothing that becomes contaminated. Wash promptly with soap & water if skin becomes contaminated.

9. PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE Colourless Liquid

ODOUR None

SOLUBILITY / MISCIBLE WITH Water; Ethanol

10. STABILITY AND REACTIVITY

STABILITY

Stable under normal temperature conditions and recommended use.

INCOMPATIBLE MATERIALS

Strong oxidizers.

11. TOXICOLOGICAL INFORMATION

INHALATION

May be harmful if inhaled. May cause respiratory tract irritation.

INGESTION

May be harmful if swallowed.

SKIN CONTACT

May be harmful if absorbed through skin. May cause skin irritation.

EYE CONTACT

May cause eye irritation.

12. ECOLOGICAL INFORMATION

ECOTOXICITY

None recorded.

Clin-Tech Ltd.

13. DISPOSAL CONSIDERATIONS

DISPOSAL METHODS

Dispose of waste and residues in accordance with local authority requirements.

14. TRANSPORT INFORMATION

	UN#	Class	Packing	Proper shipping name
RID / ADR				
<u>IATA</u>	Not regarded as hazardous for transport			
<u>IMDG</u>				

15. REGULATORY INFORMATION

STATUTORY INSTRUMENTS

This MSDS complies with regulation (EC) no. 1907/2006.

Workplace Exposure Limits EH40/2011.

16. OTHER INFORMATION

GENERAL INFORMATION

This safety data sheet covers all codes for this product.

REVISION INFORMATION

Revision No. 05

Revision Date 09/2015

Comments CHIP regulations removed

Replaces Document Generated 04/2013

HAZARD STATEMENTS IN FULL

H314 Causes severe skin burns and eye damage

H335 May cause respiratory irritation