


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (μL)	50 – 250μL	
MATERIAL	PS	
	MM	INCH
A , BASE LENGTH	127,7	5,03
B , BASE WIDTH	85,7	3,37
C , TOP LENGTH	123,3	4,85
D , TOP WIDTH	81,2	3,20
E , FLANGE LONG SIDE	2,7	0,11
F , OVERALL HEIGHT	14,5	0,57
G1 , DIA TOP	Ø6,8	Ø0,27
G2 , DIA BOTTOM	Ø6,4	Ø0,25
H , A - I LOCATION (X)	14,4	0,57
J , A - I LOCATION (Y)	11,4	0,45
K , WELL TO WELL SPACING	9	0,35
L , WELL DEPTH	10,3	0,41
M , BOTTOM OF WELL DISTANCE	4,3	0,17
N , FLANGE SHORT SIDE	2,7	0,11

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/ PK/CASE
249662	NON-TREATED	CLEAR	+	-	1/50

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



ThermoFisher
SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O.Box 280
DK-4000 Roskilde • Denmark
Tel +45 4631 2000 • Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT:

V96 MICROWELL PLATE, PS

DRAWN BY/DATE:

CHANGE ORDER NO.:

VALID FROM:

APPROVED BY/DATE:

FORM:

ITEM NUMBER:

MATERIAL:

DRAWING NO.:

VERSION:

PAGE TYPE:

SCALE:

METRIC

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC