CATALOG NUMBER		178883	
NOM. CAP. (mL)		800	
САР ТҮРЕ		FILTER	
CULTURE AREA (cm <sup>2</sup> ):		175	
MATERIAL	FLASK	PS	
	САР	BLUE HDPE	
		INCH	MM
A (NECK I.D.)		0.99	25.1
B (LENGTH W/O CAP)		7.64	194.1
C (LENGTH W/ CAP)		7.59	192.8
D (BODY WIDTH)		4.75	120.7
E (BODY DEPTH):		2.01	51.1

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.

2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

