

Eppendorf Certificate

Certificate of Purity - sterile

This package contains a high-quality consumable manufactured under the sterile Eppendorf Purity Standard.

The Eppendorf sterile consumables are produced in an ISO class 8 clean room environment according to ISO 14644-1.

For this product Eppendorf certifies the following:

- > sterile
- > pyrogen-free



Quality control and subsequent certification is done by an independent laboratory accredited to ISO 17025. Eppendorf guarantees the conformity within the following limits:

Sterility in accordance with USP, Ph. Eur. 2.6.12

Pyrogens < 0.001 EU/mL (tested according to Ph. Eur. 2.6.14 (LAL test))

(Endotoxins)

Lot-specific certificates are available on request or on the internet at www.eppendorf.com/certificates.

August, 2019

Joana Tziolis Expert Liquid Handling LH-Consumables

Your local distributor: www.eppendorf.com/contact Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany E-mail: eppendorf@eppendorf.com

ISO 9001 Certified ISO 14001 Certified

Monika Schneider Head of Global Quality Management & Regulatory Affairs

> ISO 13485 Certified