

ahn myTube MCT

PRODUCT DATA SHEET

Micro Tube 0.5 mL



Product Description

Catalogue Number	3-116-C5-0
Product Description	AHN myTube MCT, 0.5 mL, clear, sterile, bag
Sales Unit	Case / 5 x 1000 tubes

Product Characteristics

Graduations	Yes
Sterilised	Yes
Autoclavable (Max °C, Period)	Yes, 121°C, 15 min
Centrifuge Specification	Up to 20.000 x g
Temperature Range	-80°C to +121°C
Purity	Free of RNase, DNase, Human DNA, Pyrogens, Endotoxin

Size

Max. Volume	0.6 mL
Height of Product	31.6 mm

Material / Colours

Product Material	Polypropylene (PP)
Colour of Product	Clear

Packaging

Type of Packaging	Bag
Pieces / Bag	1000 tubes
Pieces / Case	5000 tubes
Dimension Case Carton (H x W x D)	399 x 273 x 400 mm

Storage and Handling

Temperature	Room Temperature, ideally 24°C
UV Light	Protect from UV light

Note: Improper storage can initiate degradation, which results in colour changes and can have negative effects on the physical properties of the products.

Safety: This product is not classified as dangerous.

