Corning[®] Stripettor[™] Ultra Pipet Controller

Designed for comfort and sample safety




This light, ergonomic, motorized pipet controller is designed for use with glass and plastic serological pipets from 0.5 to 100 mL. Conveniently positioned switches change operating modes and speeds quickly to handle different liquid volumes and various viscosities. Aspirating and dispensing speed is controlled through the finger triggers. When fully charged, the batteries enable up to 8 hours of continuous use. The large LCD display clearly indicates battery status, pipetting mode, and speed.

Features

- LCD display for continuous visualization of speed, modes, and battery status
- ▶ Colored nose pieces for user or application identification
- Well-positioned function buttons
- 3-position charging stand











Charging position Rest position

Break position





Hydrophobic PTFE filters



Battery cover with wings



Extra nose pieces



3 position charging stand

Ordering Information

Corning® Stripettor™ Ultra Pipet Controller

| Cat. No. | Description | | | | |
|-------------|--|--|--|--|--|
| 4099 | Corning Stripettor Ultra Pipet Controller, with universal power supply, three 0.2 μm filters, 3 colored nose pieces, wall hanger, and 3-position charging stand | | | | |
| Accessories | | | | | |
| 4991 | Hydrophobic PTFE filter 0.2 μm (5 pcs) | | | | |
| 4992 | Hydrophobic PTFE filter 0.45 μm (5 pcs) | | | | |
| 4993 | Silicone pipet holder | | | | |
| 4994 | Battery set, 3 x AAA, NiHM | | | | |
| 4995 | Battery cover with wings | | | | |
| 4996 | Universal power supply, 9V with set of adaptors (EU, US, UK, and AUS) | | | | |
| 4997 | Colored nose pieces, set of 4 (orange, yellow, green, red) | | | | |
| 4998 | Charging stand, 3-positions | | | | |

Related Products

Serological Pipets, Individually Wrapped, Paper/Plastic Wrap

| Cat. No. | Capacity (mL) | Graduations (mL) | Negative Graduations (mL) | Color Strip | Qty/Pk | Qty/Cs |
|----------|---------------|------------------|------------------------------|-------------|--------|--------|
| 4485 | 1 | 1/100 | 0.2 | Yellow | 50 | 1,000 |
| 4486 | 2 | 1/100 | 0.2 | Green | 50 | 1,000 |
| 4487 | 5 | 1/10 | 2.5 | Blue | 50 | 200 |
| 4488 | 10 | 1/10 | 3.0 | Orange | 50 | 200 |
| 4489 | 25 | 2/10 | 10.0 | Red | 25 | 200 |
| 4490 | 50 | 1/2 | 10.0 | Purple | 25 | 100 |
| 4491 | 100 | 1/1 | N/A | Aqua | 10 | 100 |

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit **www.corning.com/lifesciences** or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476
www.corning.com/lifesciences

ASIA/PACIFIC Australia/New Zealand t 61 427286832 Chinese Mainland t 86 21 3338 4338 f 86 21 3338 4300

India t 91 124 4604000 f 91 124 4604099 Japan t 81 3-3586 1996 f 81 3-3586 1291 Korea t 82 2-796-9500 f 82 2-796-9300 Singapore

t 65 6572-9740 f 65 6735-2913 **Taiwan** t 886 2-2716-0338 f 886 2-2516-7500

CSEurope@corning.com

EUROPE

France

All Other European Countries t +31 (0) 206 59 60 51 f +31 (0) 206 59 76 73

LATIN AMERICA grupoLA@corning.com Brazil t 55 (11) 3089-7400 Mexico t (52-81) 8158-8400